

# Collection Method of Aptamers

Application note Ref : AN 2.006-V1

Aptamers are DNA (or RNA) oligonucleotides selected from random libraries of DNA sequences. They fold into a three dimensional structure to bind different classes of targets with high affinity and selectivity. Both the chances for the aptamer to be selected and the quality of the selected aptamer are largely dependent on the method selection. In this note, we optimised the parameters of collection on a standard mixture to validate the collection method in a small range of migration times to get the slower koff.

## Capillary Electrophoresis Collection Method

In the aim to standardise a collection method, we analysed a mixture of three sugars: glucose (DP1), maltose (DP2) and maltotriose (DP3) labelled with APTS. We injected a plug of 100nL at  $c_i=10^{-7}$  mol/L and we measured each migration time. Then, we set up a collection method and we performed 8 collections. Sugars were collected separately in three vials containing 15µL of buffer (Fig.1). The collected samples were injected again to estimate the yield and the reproducibility of our collections (Fig. 2).




Figure 1: Analysis by CE-UV-LEDIF of three sugars (tm: migration time; A:area)  
 DP1  $t_m=16.4min - A1=24.9$ ; DP2  $t_m=18.7min - A2= 30.6$ ; DP3  $t_m=21.1min - A3= 29.3$

**Instruments:**

Capillary Electrophoresis: Agilent Technologies CE7100  
 Detector: Picometrics ZETALIF LED 480nm/30nm


**Sample:**

Glucose (DP1), maltose (DP2) and maltotriose (DP3) are labelled with APTS and diluted in water. A plug of 100nL is injected into the capillary.

**Method:**

- Capillary: 90cm x 50µm ID (effective length 69cm)
- Migration buffer: 40mM aminocaproic acid pH=4.5 + 0,02% hydroxypropylmethylcellulose.
- Voltage: +20kV
- Injection: 30mbars, 10s
- Cassette temperature : 25°C

We checked if the peak area of the collected products corresponds to the dilution factor DF:  
 $c_i=10^{-7}$  mol/L;  $v_i=100nL$  8 collects and  $v_f=15µL$  so  $DF = 8 \cdot 10^{-7} / 15 \cdot 10^{-6} = 0.05$ .



DP1	Glucose	→ A=1.2
DP2	Maltose	→ A= 1.8
DP3	Maltotriose	→ A=1.4

$AF1_{theoretical}=DF \cdot A11= 24.9 \cdot 0,05=1.3$	↔	$A1_{experimental}= 1.2$
$AF2_{theoretical}=DF \cdot A12= 30.6 \cdot 0,05=1.5$	↔	$A2_{experimental}= 1.8$
$AF3_{theoretical}=DF \cdot A13= 29.3 \cdot 0,05=1.5$	↔	$A3_{experimental}= 1.4$

Migration times are different due to the fact that samples were collected in a buffer.

Figure 2: Injection of products collected by CE-UV-LEDIF DP1  $t_m=16.4min - A1=1.2$ ; DP2  $t_m=18.4min - A2=1.8$ ; DP3  $t_m= 21.5min - A3=1.4$

**Conclusion:**

We are able to collect each of three sugars with a yield close to 100%. We developed an efficient method of complex separation which leads us to collect aptamers with an excellent yield and to overpass these bottleneck steps.