

# BIABooster CE-LIF SYSTEM



## **BIABooster Platform:**

- Agilent Technologies 7100 Capillary Electrophoresis
- Picometrics Zetalif™ Laser Induced Fluorescence Detector
- Capillary Device

This note describes how the  $\mu$ LAS technology is used with the Agilent Technologies 7100 Capillary Electrophoresis coupled to Picometrics Laser Induced Fluorescence Detector (CE-LIF).

The BIABooster system provides high sensitivity and extended range of DNA size analysis. It is also possible to perform all modes of capillary electrophoresis that are available with the Agilent Technologies CE.

Learn more about **Adelis**

 [info@adelis-tech.com](mailto:info@adelis-tech.com)

[www.adelis-tech.com](http://www.adelis-tech.com)

# Instrumental Set-Up


The BIABooster solution consists of the following items:

- Agilent Technologies 7100 Capillary Electrophoresis System including all functional hardware for performing CE separation
- Picometrics Zetalif™ Laser Detector including a Detector, Optical Cell, LIF Cassette, Laser light source with the corresponding emission filter block, LIF Driver for Agilent Technologies Software.
- Proprietary capillary device for  $\mu$ LAS technology
- BIABooster Analytics software to quantify, qualify and size DNA.

## About $\mu$ LAS Technology


$\mu$ LAS technology simultaneously uses electric field and pressure in a viscoelastic fluid to analyse DNA. The capillary device is specially designed to take full benefit of the technology.

**Basic principle** : DNA is subjected to a pressure-driven viscoelastic flow in combination with a counter-electrophoresis. In these conditions, DNA undergoes a viscoelastic force oriented toward the channel walls, the amplitude of which depends on its size. Because of the parabolic velocity profile of the flow, DNA molecules are transported by the fluid at a rate which depends on their size, like in gel electrophoresis.




$\mu$ LAS basically covers three functions:

**C O N C E N T R A T E**  
On-line concentration




**To concentrate DNA "in line"**:  $\mu$ LAS enables the design of "DNA traps" within a flow, by joining two capillaries of different diameters. In these traps, molecules are confined precisely at capillary junction, allowing to stack them and/or remove them from the flow. DNA can therefore be concentrated before being separated, for unmatched analytical sensitivity.

**S E P A R A T E**  
Size separation



**To separate and quantify DNA fragments** of different sizes, which have equal electrophoretic mobilities. As such, it is a new separation technique, complementary to gel electrophoresis since molecules are separated in solution.


**P U R I F Y**  
On-line purification



**To isolate and purify DNA fragments of interest**: this functionality is obtained by combining the two functions of concentration and separation, and by adding a collecting function available with the 7100 CE. DNAs of the desired size can be selected, excluding other DNAs present in the sample.

# Proprietary Capillary Device

A modified cassette accommodates the integration of a LIF detection system and any type of  $\mu$ LAS capillary device.



On-line DNA concentration at a  $\mu$ LAS capillary junction before separation

# Specifications for cfDNA Analysis

cfDNA is a promising biomarker for non-invasive monitoring of cancer disease. Getting a DNA profile of plasmatic free DNA is difficult using existing electrophoresis systems. But it is an easy thing with the BIABooster system using the **DNA1K kit** which is used for Quality Control of circulating DNA previously purified from plasma.


Analytical specifications	Value
Sizing range	0.1-1.5 kb (DNA 1K kit)
Limit of detection (S/N = 3)	10 pg/ml at 1kb – 100 pg/ml at 100 bp
Sizing accuracy*	+/- 3%
Sizing reproducibility	1% CV
Quantitative range*	5-1000 ng/ml for cfDNA
Quantitative precision	10% CV
Quantitative accuracy	20%
Minimum sample volume	10 $\mu$ l (1 $\mu$ l injected)
Maximum salt concentration	150 mM

\* Determined using a commercial 100 bp ladder as a sample, different from the standard used for reference.  
Excitation wavelength : 488nm Laser or LED option available


# Some BIABooster Applications

## cfDNA Analysis in Plasma


- Unrivalled Sensitivity
- Unrivalled Robustness
- Purification + Concentration + Separation in 30 minutes



Concentration and separation of 100 bp ladder with 100 mM NaCl (blue), 150 mM NaCl (red), 200mM NaCl (green), 250 mM NaCl (pink)



Analysis of purified cfDNA (blue) and 100 bp ladder (red)



Analysis of endogenous DNA of plasma (with Proteinase K + SDS treatment)


The BIABooster system is able to provide a profile of circulating DNA directly from **plasma**.

## cfDNA Sample Quality Control




Qualify your samples before PCR and sequencing analysis:


- Concentration
- Integrity
- Genomic DNA contamination




cfDNA at 210 ng/ml  
Typical profile



cfDNA at 270 ng/ml  
HMW contamination





cfDNA at 4 ng/ml  
Low concentration of cfDNA  
with good integrity



cfDNA at 7 ng/ml  
Low concentration of  
degraded cfDNA


# High Molecular Weight DNA Separation

Next generation long read sequencing has increased interest in large DNA fragment analysis. BIABooster platform enables to quantify and qualify DNA fragments up to 150 kb.



## cfDNA Fractionation

$\mu$ LAS technology can be used to select a DNA size range of interest. This has been used for cell free circulating DNA in which tumoral cfDNA has been reported to have a smaller size compared to constitutional cfDNA. The isolation of tumoral cfDNA is expected to provide a better sensitivity for mutation detection.



## References:

CfDNA Biomarker research applications, see [article in Anal. Chem., 2018, 90 \(6\), pp 3766–3774](#)

---

The Adelis logo is displayed in a stylized, glowing yellow font within a dark rectangular box. A thin black horizontal line extends from the left edge of the page to the left side of this box.

Learn more about **Adelis**

 [\*\*info@adelis-tech.com\*\*](mailto:info@adelis-tech.com)

[\*\*www.adelis-tech.com\*\*](http://www.adelis-tech.com)

---

Product Note No: BIA001-V2

Specifications subject to change without notice as part of our ongoing quality improvement program.  
July 2019