


## Introduction – the SChISM study


Immune-checkpoint inhibitors (ICIs) have revolutionized cancer treatment, but a large proportion of patients is resistant to the treatment. The SChISM study (SI 20.07.23.62212) was designed to monitor and predict response by measuring the size distribution of cfDNA in patients ICI-treated in monotherapy or combination.

### CLINICAL CHALLENGE

- 20-40%<sup>1</sup> of advanced cancer patients show long-term response to ICIs
- Many experience early progression (EP), progression at 1<sup>st</sup> radiological assessment
- Need for reliable early predictive biomarkers




## BIABooster System




BIABooster system is based on a dual hydro- and electro-kinetic actuation, which pushes DNA towards walls, where flow rate is reduced. It gives an electrophoresis-like view of the sample, with a ~1,000-fold increase in sensitivity over state-of-the-art CGE, thanks to in-line DNA concentration. It can also operate without DNA extraction from plasma<sup>2,3</sup>.


### Fluorescence traces reveal large fragments (healthy donors)



### Easy workflow







### Size distribution descriptors




## Results

Patients clustered according to their fragment size distribution<sup>5</sup>



### An intriguing ~3.1 kb fragment is often found in size profiles



cfDNA was extracted and analysed for fragments greater than 1.5 kb:

Left: ID-Xtract extraction (ID-Solutions) – SChISM Patients  
Right: QIASymphony extraction (Qiagen) – Healthy Donors from DIAL study (SI 22.01274.000075)

Presence of the ~3.1 kb fragment:

- Healthy donors: 11 / 13 samples
- Cancer patients: 10/27 samples – no link with EP
- Working hypotheses:
  - Artefact caused by the fixative included in the blood collection tubes (Roche BCT)
  - If not, NETOSIS or mitochondrial derived cfDNA




For more information contact: [fginot@adelis-tech.com](mailto:fginot@adelis-tech.com)

### Presence of long fragments predict no early progression<sup>4,5</sup>


Logistic and Cox regressions results: high-molecular-weight fragments ratio as the best predictor of no EP and long PFS.

(AUC: Area Under the ROC curve; HR: Hazard Ratio; NLR: Neutrophil-to-lymphocyte ratio; MV: Multivariate (adjustment with confounding features: age, sex, ECOG, NLR, pathology); OR: Odds Ratio; UV: univariate)  
(Significance: \*\*\*\*: p-value < 0.0001; \*\*\*: < 0.001; \*\*: < 0.01; \*: < 0.05)

VARIABLE	EARLY PROGRESSION		PROGRESSION-FREE SURVIVAL			
	AUC	OR UV	OR MV	C INDEX	HR UV	HR MV
CFDNA METRICS						
R <sub>5-1650</sub>	0.73	0.39 (0.25 – 0.62) ****	0.3 (0.14 – 0.63) **	0.69	0.54 (0.43 – 0.68) ****	0.43 (0.32 – 0.58) ****
C <sub>160-220</sub>	0.70	2.1 (1.4 – 3.2) **	2.2 (1.2 – 4) **	0.66	1.7 (1.4 – 2.1) ****	1.9 (1.5 – 2.5) ****
C <sub>580-1650</sub>	0.69	0.54 (0.36 – 0.8) **	0.46 (0.25 – 0.83) *	0.66	0.6 (0.48 – 0.75) ****	0.49 (0.37 – 0.64) ****
C <sub>280-320</sub>	0.67	1.8 (1.2 – 2.7) **	2.3 (1.2 – 4.3) **	0.64	1.6 (1.3 – 2) ****	1.8 (1.4 – 2.4) ****
R <sub>580-1650</sub>	0.65	0.51 (0.32 – 0.82) **	0.53 (0.25 – 1.1)	0.64	0.55 (0.4 – 0.75) ****	0.51 (0.34 – 0.75) ***
P <sub>2</sub>	0.64	0.56 (0.36 – 0.86) **	0.39 (0.19 – 0.79) **	0.63	0.61 (0.47 – 0.79) ***	0.54 (0.39 – 0.73) ****
HP <sub>1</sub>	0.66	1.7 (1.1 – 2.5) **	2 (1.1 – 3.6) *	0.63	1.5 (1.2 – 1.9) ****	1.9 (1.4 – 2.5) ****
C <sub>700</sub>	0.65	1.7 (1.1 – 2.5) **	2.2 (1.3 – 3.9) *	0.62	1.5 (1.2 – 1.9) ****	1.9 (1.4 – 2.5) ****
CLINICAL AND BIOLOGICAL VARIABLES						
NLR	0.70	2.5 (1.5 – 4.2) ***		0.63	1.6 (1.2 – 2.2) **	



Higher ratio of long cfDNA fragments is associated with extended PFS and no EP under ICIs



## Conclusions

- ICI-treated patients with lower cfDNA fragmentation tend to not early progress and to have longer PFS.
- Our best baseline cfDNA biomarker is a strong predictor and outperforms the Tumor Mutational Burden (AUC = 0.69, PPV = 0.42) and PD-L1 immunohistochemistry (AUC = 0.65, PPV = 0.34) in the prediction anti-PD-1/PD-L1 response.
- An intriguing 3.1 kb fragment is often found, which needs to be identified.
- The BIABooster system enables studying full cfDNA size distribution in a straightforward manner.

## References

1. Sharma P, et al. Primary, Adaptive, and Acquired Resistance to Cancer Immunotherapy. Cell. 2017;168(4):707-723.
2. Andriamanampisoa CL, et al. BIABooster: Online DNA Concentration and Size Profiling with a Limit of Detection of 10 fg/µL and Application to High-Sensitivity Characterization of Circulating Cell-Free DNA. Anal Chem. 2018 Mar 20;90(6):3766-3774.
3. Boutonnet A, et al. Size and concentration of cell-free DNA measured directly from blood plasma, without prior DNA extraction. Anal Chem. 2023;95(24):9263-9270.
4. Nguyen Phuong L, Salas S, Benzekry S. Computational Modeling for Circulating Cell-Free DNA in Clinical Oncology. JCO Clin Cancer Inform. 2025;(9):e2400224.
5. Nguyen Phuong L, PhD thesis, Mechanistic modeling of circulating cell-free DNA for prediction of response to immunotherapy, Aix-Marseille University, 2025.